

ImpulseAdjusting
TECHNIQUE[®]

mechanical
INNOVATIONS


ImpulseAdjusting TECHNIQUE[®]

Lower Extremity

Hip, Knee, Ankle and Foot

ImpulseAdjusting
TECHNIQUE[®]

1



Mechanical vs. Functional Instability

ImpulseAdjusting
TECHNIQUE[®]

mechanical
INNOVATIONS

Injury and various diseases may lead to two types of dysfunction: 1) mechanical instability and 2) functional instability.”


1) *Mechanical instability* is made up of the “hardware” components – muscles, joints, tendons, and ligaments. This would include conditions such as ruptured tendons and ligaments, damaged bones, torn menisci, DJD, RA, etc., and diseases causing muscle weakness (polio, CVA, etc). Many of these conditions require corrective surgery and/or are permanent. Some, such as weakness from a CVA are not amenable to surgical correction.

“Tale of the tape,” Jardine, *Advance for Physical Therapy and Rehab Medicine*, Jan 18, 2010

ImpulseAdjusting
TECHNIQUE[®]

2

1




Mechanical vs. Functional Instability

2) **Functional instability** is made up of the “software” components – the operating system or neural programming behind muscle activation and coordinated action between muscles to produce smooth movement and provide stabilization of the joint.

Conditions causing mechanical instability will always cause some functional instability.

Improvement of the functional instability will lessen the mechanical instability and improve function, but not to the degree that it was before the condition that caused the damage.

“Tale of the tape,” Jardine, *Advance for Physical Therapy and Rehab Medicine*, Jan 18, 2010



3



Frequent conditions of the ankle & foot

- **Plantar Fasciitis / Heel Pain**
- **Achilles Tendonitis**
- **Inversion Ankle Sprain**
- **Eversion Ankle Sprain**
- **Instep Pain**
- **Metatarsalgia**




4

PATIENT HISTORY


Initial Hypotheses Based on Historical Findings

Patient Reports	Initial Hypothesis
Patient reports a traumatic incident resulting in either forced inversion or eversion	Possible ankle sprain ^{1,2} Possible fracture Possible peroneal nerve involvement (if mechanism of injury is inversion) ³⁻⁵
Patient reports trauma to ankle that included tibial rotation on a planted foot	Possible syndesmotic sprain ¹
Patient notes tenderness of anterior shin and may exhibit excessive pronation. Symptoms may be exacerbated by repetitive weight-bearing activities	Possible medial tibial stress syndrome ⁶
Patient reports traumatic event resulting in inability to plantarflex ankle	Possible Achilles tendon rupture
Patient reports pain with stretch of calf muscles and during gait (toe push off)	Possible Achilles tendonitis ⁷ Possible Sever's disease ¹
Patient reports pain at heel with first few steps out of bed after prolonged periods of walking	Possible plantar fasciitis
Patient reports pain or paresthesias in plantar surface of foot	Possible tarsal tunnel syndrome ¹ Possible sciatica Possible lumbar radiculopathy
Patient reports pain on plantar surface of foot between 3rd and 4th metatarsals. Might also state that pain is worse when walking with shoes compared with barefoot	Possible Morton's neuroma ⁷ Possible metatarsalgia


5



6




Plantar Fasciitis




ImpulseAdjusting
TECHNIQUE™

7




Calcaneal Apophysitis (Seer's Disease)

Most common cause of heel pain in children and teenagers:



ImpulseAdjusting
TECHNIQUE™

8



9

Plantar Fasciitis

ImpulseAdjusting TECHNIQUE


mechanical INNOVATIONS

Combination Adjustment


1. Medial Calcaneus
2. Inferior-medial Talus
3. Inferior-medial Navicular
4. Inferior Medial Cuneiform
5. Inferior 1st Metatarsal
6. Inferior Lateral Cuneiform

ImpulseAdjusting TECHNIQUE


10




11




12




13




14




15




16



17




18



Achilles Tendon Injury

Must rule out complete rupture

Thompson Test




- 1 The patient lies in a supine position.
- 2 The examiner applies a squeeze to the calf of the patient's affected leg.
- 3 A positive test is a nonresponse during the squeeze test.


UTILITY SCORE ?

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Thompson & Doherty ²⁹	NT	40	NT	NA	NA	7

Comments: The test has surprisingly low sensitivity. Concurrent patient history is essential when performing this test.




19




Achilles Tendinitis

Combination Adjustment Protocol


1. Posterior-Superior Calcaneus
2. Posterior Distal Femur
3. Posterior Proximal Tibia



20



Achilles Tendinitis





Adjustment Directive

Posterior Superior Calcaneus


Stylus: Single

LOD: Inferior & Anterior


Setting: Medium

21



Achilles Tendinitis





Adjustment Directive

Posterior Distal Femur


Stylus: Large Dual

LOD: Anterior

Setting: Medium

22



Achilles Tendinitis


Adjustment Directive

**Posterior
Proximal Tibia**

Stylus: Large Dual

LOD: Anterior

Setting: Medium




ImpulseAdjusting
TECHNIQUE

23

“Continued sports activity using a pain-monitoring model during rehabilitation in patients with achilles tendinopathy,” Silbernagel et al, *American Journal of Sports Medicine*, 2007:35


- **Rehab for both conditions involves “heavy-load eccentric” exercises**
- **Have been shown repeatedly to be highly effective.**
- **Cortisone injections have been shown to be ineffective, dangerous, and weaken the tendon making rupture more likely.**



24



25



26



Inversion Ankle Sprain

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

- Up to 40% of all athletic injuries are related to ankle sprains.
- Ankle injury rates reach as high as 50% in basketball players and 29% in soccer players.

Smith & Reischl. Treatment of ankle sprains in young athletes. The American Journal of Sports Medicine, 1986:14

ImpulseAdjusting
TECHNIQUE™

27



Inversion Ankle Sprain

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

- **Must rule out:**
 - **Dislocation of talus out of the mortise**
 - **Fractures (including distal end of the fibula)**
 - **Complete ligament ruptures**


ImpulseAdjusting
TECHNIQUE™

28



Inversion Ankle Sprain


neuromechanical
INNOVATIONS



ImpulseAdjusting
TECHNIQUE

29

Anterior Drawer Test



- 1 The patient lies in a supine position. The ankle is prepositioned into slight plantarflexion.
- 2 The examiner provides an anterior glide of the calcaneus and talus on the stabilized tibia.
- 3 A positive test is excessive translation of one side in comparison to the opposite extremity.

UTILITY SCORE **2**

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Hertel et al. ¹³	NT	78	75	3.1	0.29	8
Phisitkul et al. ²²	NT	100	100	Inf	Inf	7

Comments: The test is designed to measure damage to the anterior talofibular ligament. The examiner should observe the presence of a dimple or sulcus sign near the region of the anterior talofibular ligament. Phisitkul's study used cadavers and was poorly performed.

30

Inversion Ankle Sprain

ImpulseAdjusting
TECHNIQUE™

neuromechanical
INNOVATIONS

Talar-Tilt Sign
(Test for tear of calcaneofibular and anterior talofibular ligaments)

Examiner firmly rotates foot in varus. Tear of calcaneofibular ligament permits excessive mobility in this direction (leg firmly fixed by other hand)

Anterior talofibular ligament-torn

Calcaneofibular ligament-torn

ImpulseAdjusting
TECHNIQUE™

31

Anterior Lateral Drawer Test

- 1 The patient assumes a sitting position.
- 2 The examiner stabilizes the lower leg just above the ankle. The other hand provides an anterior directed force, measurement of talar translation, and control of ankle plantarflexion.
- 3 A positive test is 3 millimeters or more of translation.

UTILITY SCORE 3

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Phisitkul et al. ²²	NT	75	50	1.5	0.50	7

Comments: The study used cadavers and was poorly performed.

32



Figure 9. Illustration depicting the talar tilt test to assess stability of the calcaneofibular ligament. (Reprinted with permission from Hoppenfeld S, Hutton R. Physical examination of the spine and extremities. New York: Appleton-Century-Crofts; 1976:227.)

33

Medial Talar Tilt Stress Test



- 1 The patient is placed in a sitting or supine position.
- 2 The examiner grasps the ankle of the patient at the malleoli.
- 3 The examiner applies a quick medial thrust to the calcaneus.
- 4 A positive test is excessive laxity when compared to the opposite side.

UTILITY SCORE 2

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Hertel et al. ¹³	NT	67	75	2.7	0.44	8
Comments: Expect positive findings after inversion sprains.						

34



Fibular Translation Test



- 1 The patient lies in a sidelying position.
- 2 The examiner applies anterior and posterior forces on the fibula at the level of the syndesmosis.
- 3 A positive test is pain during translation and more displacement to the fibula than the compared side.

UTILITY SCORE **2**

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Beumer et al. ³	NT	82	88	6.8	0.2	8
Comments: Beumer et al. ³ only found increased translation when all ligaments were removed in cadavers.						

35



ImpulseAdjusting
TECHNIQUE™

Inversion Ankle Sprain

ImpulseAdjusting
mechanical
INNOVATIONS


Combination Adjustment

1. Anterior Lateral Superior Talus
2. Posterior Distal Fibula


ImpulseAdjusting
TECHNIQUE™

36

18



Inversion Ankle Sprain





Adjustment Directive

Anterior Lateral Talus


Stylus: Single

LOD: Posterior Medial


Setting: Medium

37



Inversion Ankle Sprain





Adjustment Directive

Anterior Lateral Talus


Stylus: Single

LOD: Posterior Medial


Setting: Medium

38





Inversion Ankle Sprain




Adjustment Directive

Posterior Distal Fibula


Stylus: Single
LOD: Anterior
Setting: Medium

39





Inversion Ankle Sprain



Adjustment Directive

Posterior Distal Fibula

Stylus: Single
LOD: Anterior
Setting: Medium

40

Inversion Ankle Sprain Rehab

Alphabet Ankle



Seated with leg out in front.
Place towel under lower leg and relax ankle.
Move ankle so as to draw the letters of the alphabet A,B,C...
Try drawing both upright figure 8's and side lying 'lazy 8's in both directions.



41

Inversion Ankle Sprain Rehab

Step Down



Standing on a low step. Step down with the unaffected leg. Raise body weight back up onto step and repeat.



42

Exercises

From: norespond@3drx.com
Sent: Thu 10/04/12 9:09 AM
To: docjim47@live.com



Chiro-Plus

Has generated a report on 10/4/2012

You can also view report in details, [Click](#)

[view report](#)

Rehab Section

Step Down - FLO206

Standing on a low step, step down slowly with one leg. Control the lowering with upper leg. Upper leg then lifts body weight back up, foot returns to step. Repeat. Maintain an abdominal hollow. Stop if the exercise causes pain.



43

Proprioceptive Training

A-P Rock



Keeping finger tips on wall or chair for support, step onto rocker board with one foot.

Keeping eyes looking straight forward, raise opposite foot onto board.


Maintain a short foot and readiness position while on the board.

With control, slowly rock board forward and back.




44

22



45



Eversion Ankle Sprain

mechanical INNOVATIONS

Combination Adjustment

1. Inferior Medial Talus
2. Inferior Medial Navicular
3. Inferior Medial Tibia

ImpulseAdjusting
TECHNIQUE

46

Eversion Ankle Sprain

Common Associated Ankle/Leg Subluxations

- Inferior-medial Talus
- Inferior-medial Navicular
- Inferior medial Tibia



47




Inferior-medial Talus

- LOD: Superior & Lateral
- Setting: Low or Medium


Eversion Ankle Sprain

48


24



49




50



Ankle Impingement

Forced Dorsiflexion Test




- 1 The patient assumes a sitting position.
- 2 The examiner stabilizes the distal aspect of the tibia and places his or her thumb on the anterolateral aspect of the talus near the lateral gutter. Pressure is applied.
- 3 The examiner applies a forceful dorsiflexion movement.
- 4 A positive test is reproduction of pain at the anterolateral aspect of the foot during forced dorsiflexion.


UTILITY SCORE 2

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Alonso et al. ¹	0.36 kappa	NT	NT	NA	NA	NA
Molloy et al. ¹⁷	NT	95	88	7.9	0.06	8

Comments: Alonso et al.¹ tested for a syndesmosis injury. Although the diagnostic values for the test are strong, the quality of the study and the reliability among examiners is poor.




51



Ankle Impingement

Clinical Prediction Rule of Impingement




Five of six symptoms below are considered positive for anterior ankle impingement:

- 1 Anterolateral ankle joint tenderness.
- 2 Anterolateral ankle joint swelling.
- 3 Pain with forced dorsiflexion.
- 4 Pain with single-leg squat on the affected side.
- 5 Pain with activities.
- 6 Absence of ankle instability.


UTILITY SCORE 2

Study	Reliability	Sensitivity	Specificity	LR+	LR-	QUADAS Score (0-14)
Liu et al. ¹⁵	NT	94	75	3.8	0.08	7


Comments: Some disagreement exists whether absence of ankle instability should be a rule for impingement. The quality of the single study is suspect.



52




Ankle Impingement




Combination Adjustment

1. Anterior-Superior Talus
2. Superior Navicular
3. Superior Medial Cuneiform
4. Anterior Distal Tibia




53



Anterior-Superior Talus
 -LOD: Posterior-inferior
 -Setting: Low to Medium

Instep Pain


Superior Navicular
 LOD: inferior



54



55



56



57



58

 **Anterior Tibia** 


LOD: posterior



59

Foot/Toe Conditions

- Hallux Valgus/1st MTP Sprain
- Metatarsalgia

Experience the Next Generation in Chiropractic™ **NEUROMECHANICAL INNOVATIONS®** 


60



61



62




63



64



65



66


 **Inferior MTP**
LOD: superior
Setting: Low or medium

mechanical INNOVATIONS




67


Cuboid sprain:
LOD: inferior




68



69



70




Knee

mechanical INNOVATIONS


Knee Articulations

1. Tibiofemoral Joint
2. Tibiofibular Joint
3. Patellofemoral Joint



ImpulseAdjusting
TECHNIQUE

71




Knee – Tibiofemoral Joint

mechanical INNOVATIONS

Physical Examination

1. Orthopedic Examination
2. Postural Assessment



ImpulseAdjusting
TECHNIQUE

72



Knee

mechanical INNOVATIONS


Knee Adjustment Considerations

1. Anterior
2. Tibiofibular Joint
3. Patellofemoral Joint

Impulse Adjusting TECHNIQUE

73

Occurrence of anterior cruciate ligament (ACL) injury expressed as (A) a percentage of all injuries and (B) the rate per 1000 exposures (games and practices combined, 1988-9 through 2003-4).



P Renstrom et al. Br J Sports Med 2008;42:394-412

BJSM

Copyright © BMJ Publishing Group Ltd & British Association of Sport and Exercise Medicine. All rights reserved.

74

37

The Missing ACL Link

Getting a handle on the most dangerous trend facing today's athlete **By Steve Grosserode, DPT, OCS, and Jared Vagy, DPT, OCS**

It is estimated that there are over 250,000 anterior cruciate ligament (ACL) injuries per year in the United States. More than 100,000 ACL reconstructions are performed annually, costing up to \$17 billion each year.¹

ACL tears have become an epidemic. There are more young athletes suffering from this injury than ever before.² Over half of ACL tears occur in athletes 15-20 years old.³

These injuries have serious effects on the lives of young athletes. They can lead to a decreased grade point average, depression, fear about returning to their sport, and even the loss of a collegiate scholarship. Long-term risks include up to a 70% chance of early osteoarthritis 10-15 years following injury.⁴⁻⁷

Fortunately, over the last decade, an abundance of evidence has shown that ACL injuries can be prevented with a focus on proper movement technique and proximal hip strength. Mandelbaum et al. showed that up to 88% of ACL injuries can be prevented with education, stretching, strengthening, plyometrics and sports-specific agility drills with an emphasis on proper technique, including video feedback movement training.⁸

With ACL injuries occurring at alarming rates, it is essential to develop strategies to reduce them. Athletes who are predisposed to ACL injury must be identified through screenings that analyze movement patterns and test the strength of the muscles that prevent injuries. Following the screening, a physical therapist must devise an individual training program that combines the findings from the screening with the latest research.

Poor Movement and Weak Hips
The ACL is most stressed when the knee is placed in a valgus position and the femur is internally rotated on a fixed tibia.^{9,11} The knee is placed in these harmful positions during planting, cutting and turning movements that are common in most sports. The two main factors that protect the ACL from valgus stress are proper movement patterns and muscle training — specifically proximal hip strengthening.

Movement patterns originate in the brain and are carried out via motor control. Motor control occurs when the brain and nervous systems coordinate with muscles to perform a task. Therefore, movement retraining is separate from muscle training.

The gluteus maximus and gluteus medius are the primary muscles in the hip that control abduction and internal rotation of the femur. Weakness of these muscles can lead to knee

16 ADVANCE for Physical Therapy & Rehab Medicine | December 8, 2014 | www.advanceweb.com/pt

75

Drop vertical jump test.

P Renstrom et al. Br J Sports Med 2008;42:394-412

Copyright © BMJ Publishing Group Ltd & British Association of Sport and Exercise Medicine. All rights reserved.


BJSM

76

38


Inappropriate neuromuscular control during landing:

Barber-Westin et al, *Journal of Strength and Conditioning Research*, Nov, 2010:24




77


“An integrated approach to change the outcome. Part II: Targeted neuromuscular training techniques to reduce identified ACL injury risk factors,” Myer et al, *Journal of Strength and Conditioning Research*, Aug, 2012




78




79




80



81



82




83

Inappropriate neuromuscular control during landing:
 Barber-Westin et al, *Journal of Strength and Conditioning Research*, Nov, 2010:24


a

b

84



85



86



Medial Knee Pain

mechanical
INNOVATIONS

Adjustment Directive

Medial Meniscus

Stylus: Single

LOD: Lateral

Setting: 2 (Medium)



ImpulseAdjusting
TECHNIQUE

87



Lateral Knee Pain

mechanical
INNOVATIONS

Adjustment Directive

Lateral Meniscus

Stylus: Single


LOD: Medial

Setting: 2 (Medium)




ImpulseAdjusting
TECHNIQUE

88




Posterior Lateral Knee Pain



Adjustment Directive


Posterior Proximal Fibula

Stylus: Single
LOD: Anterior
Setting: Medium




ImpulseAdjusting
TECHNIQUE

89




Posterior Lateral Knee Pain



Adjustment Directive

Posterior Proximal Fibula

Stylus: Single
LOD: Anterior
Setting: Medium



ImpulseAdjusting
TECHNIQUE

90

Patellofemoral Pain Syndrome (PFPS) (a.k.a. – Chondromalacia of the patella)

- “The major mechanism responsible for PFPS is probably forces causing a lateralization of the patella.”

“Upward squatting in individuals with and without PFPS: a biomechanical study,” Dionisio et al, *Journal of Strength and Condition Research*, May, 2011

- “Patellar malalignment has always been thought to be a significant cause of anterior knee pain.”

Hammer, “Treatment for Patellofemoral Pain Syndrome,” *Dynamic Chiropractic*, Feb 12, 2001


- A tight lateral retinaculum tilts the patella, leading to increased pressure on the lateral facet and reduces normal lateral to medial.

Ficat, et al, “Chondromalacia patellae: A system of classification,” *Clinical Orthopedics*, 144:1979

91

-SCP: Inferior Patella
-LOD: superior
-Setting: Medium or Low

Loss of superior
patellar glide



impulse
ADJUSTING INSTRUMENT

92

46

-SCP: Inferior Patella
-LOD: superior
-Setting: Medium or Low


Loss of superior patellar glide




93

-SCP: Superior Patella
-LOD: Inferior
-Setting: Medium or Low


Superior patellar pain



94



95




96

- “The pain from PFPS causes reflex inhibition of the quadriceps leading to decreased activation and atrophy.”
- The greater the degree of knee flexion (during the descent portion of the squat before the upward segment of the squat) the greater the torque and compression forces against the patellofemoral joint leading to pain and reflex inhibition.
- Thus initial strengthening of the quads should be performed over the final 40 degrees of knee extension during the upward squat.

“Upward squatting in individuals with and without PFPS: a biomechanical study,” Dionisio et al, *Journal of Strength and Condition Research*, May, 2011

97



98

49

“Comparison of manual therapy and exercise therapy in OA of the hip,” Hoeksma et al, *Arthritis and Rheumatism*, 2004:51


- Manual therapy was shown to be superior to exercise.

99


“Effect of therapeutic exercise for hip OA pain,” Hernandez-Molina et al, *Arthritis and Rheumatism*, 2008:59

- Therapeutic exercise, especially with an element of strengthening, is an efficacious treatment for hip OA.


100




101




102




The Hip Joint: Kinesiology




Most of the muscles acting on the hip function as a 3rd class levers (where the force is between the axis of rotation and the resistance) favoring large ROM and speed.




Abduction




103




The Hip Joint: Kinesiology



- At the hip there are 7 two-joint muscles that have one action at the hip and another at the knee.
- The body part that moves the most (hip or knee) will be the part least stabilized.



104



105

The iliopsoas

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS


In most individuals, the hip flexors remain strong and rarely need strengthening exercises but the hip extensors (and abductors and abdominals) become weak.

However, the hip flexors often become shortened due to excessive sitting (hips flexed positioned). Muscles habitually placed in relaxed (shortened) postures will become shorter due to the action of the muscle spindles


(spindles are facilitory mechanoreceptors that react to changes in muscle length to protect the muscle and will “take up the slack” to shorten a muscle that is habitually relaxed in order to maintain the muscle “tone.”)

ImpulseAdjusting
TECHNIQUE™


106




107




108




109



110




111



Gluteus Maximus

ImpulseAdjusting
TECHNIQUE™
mechanical
INNOVATIONS

- **During ordinary walking the Gluteus Maximus comes into action when the femur exceeds 15 degrees of extension, thus it is not used extensively in walking. (The posterior fibers of the G medius and the hamstrings are more active for hip extension during walking).**
- **It is also a powerful external rotator of the hip.**
- **It is a powerful stabilizer of the SI joint (pulls the sacrum laterally to close the SI joint)**
- **It is strengthened by applying resistance to it starting in a hip flexion position and moving into extension with the knee in flexion (the return from the squat position).**



112

The Effect of Sacroiliac Joint Pain on Muscle Recruitment

Posted on September 8, 2014 by admin

Like 321

Share

Send

Tweet 8

+1 4

Share 15

tumblr

Email

Share

Research Review: Altered Lumbopelvic Hip Muscle Recruitment in Individuals with Sacroiliac Joint Pain

By **Stefanie DiCarrado** DPT, PT, NASM CPT & CES

Edited by Brent Brookbush DPT, PT, MS, PES, CES, CSCS, ACSM H/FS

Original Citation: Hungerford, B., Gilleard, W., Hodges, P. (2003) Evidence of altered lumbopelvic muscle recruitment in the presence of sacroiliac joint pain. *Spine 28(14)*, 1593-1600 - ARTICLE

113

Results

- SIJP (symptomatic side vs CON): IO, LM, GMax onset significantly delayed in comparison to CON; IO and LM, delay greater than 20ms after initiation of movement (not anticipatory). BF onset occurred significantly earlier on symptomatic side than in CON with increased activity before and after initiation of movement
- SIJP (asymptomatic side vs CON): IO and LM significantly delayed but onset was within 20ms; no significant differences in onset of other muscles compared with controls
- SIJP (symptomatic vs asymptomatic side): Significant delays in IO, LM, and GMax onset on symptomatic side but no differences in BF, GMed, TFL, and AL.
- CON: IO and LM onset prior to initiation of movement; BF, GMed, TFL, AL, GMax onset after initiation of movement with no significant differences between R & L sides; first to fire were IO and LM followed by AL; BF decreased in activity as compared to quiet standing.

Normal feedforward response \Rightarrow

Conclusions

The presence of SIJP may alter recruitment strategies of the lumbopelvic hip stabilizers on both the symptomatic and asymptomatic sides. With continued pain and altered strategies, an individual can develop poor motor control of these muscles leading to further dysfunction, pain, and injury.

114



Sacroiliac Adjustment



Sacroiliac Joint Adjustment

SCP: Sacroiliac joint, medial aspect of PSIS

LOD: Anterior-Superior



Setting: High or Medium






115


The external rotators of the hip:


The six deep external rotators, along with the internal rotators, help maintain the correct position of the femoral head in the acetabulum.

They counter-act internal rotation produced in the hip during motions such as throwing a baseball.


116




117




118



Step 3: Hip Rotation Tests




External (posterior) rotation




Internal (anterior) rotation


Ilium on the side being examined must be held down firmly to prevent rotation of the pelvis.




119




Step 3: Hip Rotation Test



120



121



122

Hibb's test:


Step 3

Pain in low back or SI joint indicates a lumbar or SI joint lesion:

Pain in the hip indicates a hip lesion (DJD, inhibited external rotators or psoas)




123



Posterior Hip Adjustment

(For Loss of Hip Internal Rotation)






Posterior Hip Adjustment (to restore hip internal rotation)


SCP: Posterior aspect of the greater trochanter at the tendinous insertions of the 6 deep rotators, the posterior fibers of the G. medius, and the G. max.

LOD: Anterior


Setting: High or Medium





124



125



126


GLUTEUS MAXIMUS

Origin: Posterior gluteal line of the ilium and portion of the bone superior and posterior to it, posterior surface of the lower part of the sacrum, side of the coccyx, aponeurosis of the erector spinae, sacrotuberous ligament and gluteal aponeurosis.

Insertion: Larger proximal portion and superficial fibers of the distal portion of the muscle into the iliotibial tract of the fascia lata. Deep fibers of the distal portion into the gluteal tuberosity of the femur.


Action: Extends and laterally rotates the hip joint. Lower fibers assist in adduction of the hip joint; upper fibers assist in abduction. Through its insertion into the iliotibial tract, helps to stabilize the knee in extension.

Nerve: Inferior gluteal, L5, S1, 2.





Sacrum
Sacrotuberous lig.

Iliotibial band




127

Failure to normalize or greatly improve hip internal rotation restriction:

Psoas Facilitation/Inhibition:



128



Psoas syndrome:

ImpulseAdjusting
TECHNIQUE®

mechanical
INNOVATIONS

- **Psoas syndrome can produce the following:**
- Excessive anterior tilt of the pelvis;
- Anterior rotation of one ilium causing counternutation of the SI complex;
- Reduced internal rotation of the femur;
- Pubic symphysis/groin pain;
- Excessive lumbar lordosis (bilateral) or lumbar lateral curvature with vertebral rotation (unilateral);
- Inability or pain when attempting to stand erect after sitting.

ImpulseAdjusting
TECHNIQUE®


129




Psoas Test:

ImpulseAdjusting
TECHNIQUE®



mechanical
INNOVATIONS




130




Psoas Adjustment:


Contact just lateral and inferior to pubic ramus on the tendon of the psoas.



131




Psoas Test:




- **Remission of psoas syndrome most strongly predicted LBP response.**
(Restricted sacral nutation was most strongly correlated with other biomechanical dysfunctions.)

Licciardone et al, "Changes in biomechanical dysfunction and low back pain reduction with osteopathic manual treatment: Results from the Osteopathic Trial," *Manual Therapy*, April, 2014



132

The pectineus adjustment:



Contact just inferior to the pubic ramus between the ramus and the medial aspect of the superior femur

Experience the Next Generation in Chiropractic™

NEUROMECHANICAL INNOVATIONS®



133



ImpulseAdjusting
TECHNIQUE™

“The Hip Joint: Myofascial and Joint Patterns,”
Marc Heller, DC, *Dynamic Chiropractic*, May 7, 2007:

NEUROMECHANICAL
INNOVATIONS

- Internal rotation is the first motion affected with hip capsular problems and DJD.
- “When you find a severe lack of internal rotation in the hip and this motion cannot be restored, the patient is usually on their way to hip replacement surgery.”

ImpulseAdjusting
TECHNIQUE™

134

67

Reduced external rotation of right hip:

Step 3

Hip Rotation Tests

External Rotation


50-60 degrees

Positive

Hyperactive or Inhibited Internal Rotators

Adjust Anterior Trochanter

135



136



Anterior Hip Adjustment (For Loss of Hip External Rotation)


Anterior Hip Adjustment (to restore hip external rotation)

SCP: Anterior aspect of the greater trochanter at the tendinous insertions of the TFL, anterior fibers of the G. medius, and the G. minimus.


LOD: Posterior

Setting: High or Medium









137




Hip Internal Rotators: TFL

Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins








138




Hip Internal Rotators
Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins





Gluteal aponeurosis


GLUTEUS MEDIUS
Anterior fibers




GLUTEUS MINIMUS




139



Unilateral Hip Rotation Range of Motion Asymmetry in Patients With Sacroiliac Joint Regional Pain
Cibulka et al. Spine: Volume 23(9) 1 May 1998 pp 1009-15



“The presence of such asymmetry in patients with low back pain may help identify those with sacroiliac joint dysfunction.”



140