


NEUROMECHANICAL INNOVATIONS
 EXPERIENCE THE NEXT GENERATION OF CHIROPRACTIC®



Orthopedic Physical Examination Tests
An Evidence-Based Approach
 Second Edition





Chad E. Cook | Eric J. Hegedus

Introduction to Diagnostic Accuracy

Physical Examination:


- 1) Initial patient screening
- 2) Diagnostic test
- 3) Outcome measure

*Must be reliable and valid

neuromechanical.com


ImpulseAdjusting
 SYSTEM®


1


NEUROMECHANICAL INNOVATIONS
 EXPERIENCE THE NEXT GENERATION OF CHIROPRACTIC®



The Pelvis Region

neuromechanical.com


ImpulseAdjusting
 SYSTEM®

2



Physical Examination


To demonstrate a subluxation based on the physical examination, two of the following four PART criteria are required, one of which must be either Asymmetry/misalignment or Range-of-motion abnormality:

- Pain/tenderness (P)
- Asymmetry/misalignment (A)
- Range-of-motion abnormality (R)
- Tissue, tone changes (T)

3



4



Symptoms of a Sacroiliac Subluxation/Fixation:

- Low back pain/stiffness
- Hip pain/stiffness
- Groin pain
- Knee, leg, ankle pain/stiffness
- “Catch” in the hip
- Knee wants to “give out”
- Ankle “turns easy”
- Mimic disc herniation with radicular symptoms

Schamberger, *The Malalignment Syndrome: Implications for Medicine & Sport*, 2002, Elsevier
Bernard & Kirkaldy-Willis, *Clinical Orthopaedics and Related Research*, 1987: 217(266-280

ImpulseAdjusting
TECHNIQUE

5

Back Index

Chicago, Form 01-100

Patient Name _____ Date _____

This questionnaire will give your provider information about how your back condition affects your everyday life. Please answer every section by marking the one statement that applies to you. If two or more statements in one section apply, please mark the one statement that most closely describes your problem.

Oswestry Disability Questionnaire:

0-20 = mild
21-40 = moderate
41 & up = severe

| | |
|--|--|
| <h4>Pain Intensity</h4> <p><input type="radio"/> The pain comes and goes and is very mild.</p> <p><input type="radio"/> The pain is mild and does not vary much.</p> <p><input type="radio"/> The pain comes and goes and is moderate.</p> <p><input type="radio"/> The pain is moderate and does not vary much.</p> <p><input type="radio"/> The pain comes and goes and is very severe.</p> <p><input type="radio"/> The pain is very severe and does not vary much.</p> | <h4>Personal Care</h4> <p><input type="radio"/> I do not have to change my way of washing or dressing in order to avoid pain.</p> <p><input type="radio"/> I do not normally change my way of washing or dressing even though it causes some pain.</p> <p><input type="radio"/> Washing and dressing increases the pain but I manage not to change my way of doing it.</p> <p><input type="radio"/> Washing and dressing increases the pain and I find it necessary to change my way of doing it.</p> <p><input type="radio"/> Because of the pain I am unable to do some washing and dressing without help.</p> <p><input type="radio"/> Because of the pain I am unable to do any washing and dressing without help.</p> |
| <h4>Sleeping</h4> <p><input type="radio"/> I get no pain in bed.</p> <p><input type="radio"/> I get pain in bed but it does not prevent me from sleeping well.</p> <p><input type="radio"/> Because of pain my normal sleep is reduced by less than 25%.</p> <p><input type="radio"/> Because of pain my normal sleep is reduced by less than 50%.</p> <p><input type="radio"/> Because of pain my normal sleep is reduced by less than 75%.</p> <p><input type="radio"/> Pain prevents me from sleeping at all.</p> | <h4>Lifting</h4> <p><input type="radio"/> I can lift heavy weights without extra pain.</p> <p><input type="radio"/> I can lift heavy weights but it causes extra pain.</p> <p><input type="radio"/> Pain prevents me from lifting heavy weights off the floor.</p> <p><input type="radio"/> Pain prevents me from lifting heavy weights off the floor, but I can manage if they are conveniently positioned (e.g., on a table).</p> <p><input type="radio"/> Pain prevents me from lifting heavy weights off the floor, but I can manage light to medium weights if they are conveniently positioned.</p> <p><input type="radio"/> I can only lift very light weights.</p> |
| <h4>Sitting</h4> <p><input type="radio"/> I can sit in any chair as long as I like.</p> <p><input type="radio"/> I can only sit in my favorite chair as long as I like.</p> <p><input type="radio"/> Pain prevents me from sitting more than 1 hour.</p> <p><input type="radio"/> Pain prevents me from sitting more than 1/2 hour.</p> <p><input type="radio"/> Pain prevents me from sitting more than 10 minutes.</p> <p><input type="radio"/> I avoid sitting because it increases pain immediately.</p> | <h4>Traveling</h4> <p><input type="radio"/> I get no pain while traveling.</p> <p><input type="radio"/> I get some pain while traveling but none of my usual forms of travel make it worse.</p> <p><input type="radio"/> I get extra pain while traveling but it does not cause me to seek alternate forms of travel.</p> <p><input type="radio"/> I get extra pain while traveling which causes me to seek alternate forms of travel.</p> <p><input type="radio"/> Pain restricts all forms of travel except that done while lying down.</p> <p><input type="radio"/> Pain restricts all forms of travel.</p> |
| <h4>Standing</h4> <p><input type="radio"/> I can stand as long as I want without pain.</p> <p><input type="radio"/> I have some pain while standing but it does not increase with time.</p> <p><input type="radio"/> I cannot stand for longer than 1 hour without increasing pain.</p> <p><input type="radio"/> I cannot stand for longer than 1/2 hour without increasing pain.</p> <p><input type="radio"/> I cannot stand for longer than 10 minutes without increasing pain.</p> <p><input type="radio"/> I avoid standing because it increases pain immediately.</p> | <h4>Social Life</h4> <p><input type="radio"/> My social life is normal and gives me no extra pain.</p> <p><input type="radio"/> My social life is normal but increases the degree of pain.</p> <p><input type="radio"/> Pain has no significant effect on my social life apart from limiting my more energetic interests (e.g., dancing, etc).</p> <p><input type="radio"/> Pain has restricted my social life and I do not go out very often.</p> <p><input type="radio"/> Pain has restricted my social life to my home.</p> <p><input type="radio"/> I have hardly any social life because of the pain.</p> |
| <h4>Walking</h4> <p><input type="radio"/> I have no pain while walking.</p> <p><input type="radio"/> I have some pain while walking but it doesn't increase with distance.</p> <p><input type="radio"/> I cannot walk more than 1 mile without increasing pain.</p> <p><input type="radio"/> I cannot walk more than 1/2 mile without increasing pain.</p> <p><input type="radio"/> I cannot walk more than 1/4 mile without increasing pain.</p> <p><input type="radio"/> I cannot walk at all without increasing pain.</p> | <h4>Changing degree of pain</h4> <p><input type="radio"/> My pain is rapidly getting better.</p> <p><input type="radio"/> My pain fluctuates but overall is definitely getting better.</p> <p><input type="radio"/> My pain seems to be getting better but improvement is slow.</p> <p><input type="radio"/> My pain is neither getting better or worse.</p> <p><input type="radio"/> My pain is gradually worsening.</p> <p><input type="radio"/> My pain is rapidly worsening.</p> |

Back Index Score

Signature of Patient _____

6



Dreyfuss et al, "The Value of Medical History and Physical Examination in Diagnosing Sacroiliac Joint Pain,"
Spine 1996; 21(22):2594-602.

ImpulseAdjusting
 TECHNIQUE™

mechanical
 INNOVATIONS


- 12 SIJ tests evaluated by a multidisciplinary expert panel to identify a SIJ test that is clinically valuable.
- Included: Gillet test, Thigh Thrust, Patrick's Test, Gaenslen's Test, Midline Sacral Thrust, Sacral Sulcus Tenderness, Joint Play and others.
- Results compared to SI block ("gold" standard):
- 4 tests faired best in sensitivity:
 - Sacral sulcus tenderness
 - Pain over the SIJ
 - Buttock Pain
 - Patient pointing to SIJ



7

TESTS FOR SACROILIAC PAIN ORIGIN

Thigh Thrust (also known as the Ostgaard Test, 4P Test, Sacrotuberous Stress Test, and POSH Test)



- 1 The patient is positioned in supine. Resting symptoms are assessed.
- 2 The examiner stands opposite the painful side of the patient.
- 3 The hip on the painful side is flexed to 90 degrees.
- 4 The examiner places his or her hand under the sacrum to form a stable "bridge" for the sacrum.
- 5 A downward pressure is applied through the femur to force a posterior translation of the innominate. The patient's symptoms are assessed to determine if they are concordant.
- 6 A positive test is concordant pain that is posterior to the hip or near the sacroiliac joint. A positive test requires reproduction of pain on the thrust side (the side of the loaded femur).

UTILITY SCORE **2**

ImpulseAdjusting
 SYSTEM™

8

UTILITY SCORE **2**

| Study | Reliability | Sensitivity | Specificity | LR+ | LR- | QUADAS Score (0-14) |
|---------------------------------------|-------------------------|--------------------|-------------|------|-----------|---------------------|
| Laslett & Williams ²⁵ | 0.82 | NT | NT | NA | NA | NA |
| Dreyfuss et al. ¹³ | 0.64 | 36 | 50 | 0.72 | 1.28 | 10 |
| Kokmeyer et al. ²³ | 0.67 | NT | NT | NA | NA | NA |
| Damen et al. ¹¹ | NT | 62 | 72 | 2.2 | 0.53 | 8 |
| Ostgaard & Andersson ³⁶ | NT | 80 | 81 | 4.21 | 0.25 | 5 |
| Broadhurst & Bond ⁵ | NT | 80 | 100 | NA | NA | 9 |
| Albert et al. ¹ | 0.70 | 84-93 [*] | 98 | 46.5 | 0.07-0.16 | 7 |
| Laslett et al. ²⁴ | NT | 88 | 69 | 2.8 | 0.17 | 12 |
| Arab et al. ² | 0.60 right 0.40 left | NT | NT | NA | NA | NA |
| Ozgocmen et al. ³⁷ (Right) | NT | 55 | 70 | 1.91 | 0.62 | 10 |
| Ozgocmen et al. ³⁷ (Left) | NT | 45 | 86 | 3.29 | 0.63 | 10 |
| Gutke et al. ¹⁸ | NT | 88 | 89 | 8.0 | 0.13 | 7 |

Comments: One of the few sacroiliac tests that exhibits fair sensitivity. To accurately perform the test make sure the thigh is held in neutral adduction and at 90 degrees of flexion.
^{*}Ozgocmen et al.³⁷ assessed patients with acute sacroiliitis.



7 Step IAS Pelvic Analysis



ImpulseAdjusting
TECHNIQUE


mechanical
INNOVATIONS

Consideration of PART During:



1. Postural Observation
2. Prone Hip Extension Test (Active & Passive)
3. Bilateral Knee Flexion Test
4. Unilateral Knee Flexion Test
5. Hip Rotation Tests
6. Lower Extremity Muscle Compression Tests
7. Lumbo-pelvic Muscle Compression Tests





 **Step 1. Postural Asymmetry** 



11

 **Step 2. Prone Hip Extension Test** 

a  **b** 

ImpulseAdjusting
SYSTEM[®]

12

 **Step 3. Bilateral Knee Flexion Test** 



Note: this is not a knee ROM test. Both knees are flexed to the point of rectus femoris “resistance.”




13

 **Step 3. Bilateral Knee Flexion Test** 






14



15



16




Step 4. Unilateral Knee Flexion Test


Impulse **mechanical**
TECHNIQUE INNOVATIONS

ImpulseAdjusting
TECHNIQUE

17



18



19

Location of the short posterior and long posterior Sacroiliac ligaments (superior & inferior sulci):

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS


- Patients often can identify location with one finger.
- Palpation may be the most reliable indication of SI dysfunction.
- A positive Fortin finger test is a quick and simple way of identifying an SI lesion.

The anatomical diagram shows the sacroiliac joint with the following labels: Short post. sacroiliac lig., Long post. sacroiliac lig., Sacrotuberous ligament, and Superficial post. sacrospinous ligament. Two black arrows point to the superior and inferior sulci of the sacroiliac joint.


Fortin & Falco, "The Fortin finger test," *American Journal of Orthopedics*, 1997:26(7)

ImpulseAdjusting
TECHNIQUE™

20




SI joint Vectors – coronal section



SI capsule contains a dense plexus of unmyelinated nerve fibers indicative of a nociceptive receptor system.

Wyke, "Receptor systems in lumbosacral tissues in relation to the production of low back pain," *American Academy of Orthopaedic Surgeons Symposium on Idiopathic Low Back Pain*, Mosby, 1982



21



Impulsive Sacroiliac Adjustment

Sacroiliac Joint Adjustment

SCP: Sacroiliac joint, medial aspect of PSIS



LOD: Anterior-Superior


Setting: High (or Medium in special circ.)







22



 **Step 5. Hip Rotation Tests** 





Ilium on the side being examined must be held down firmly to prevent rotation of the pelvis.



23


 **Step 5. Hip Rotation Tests** 






24


Flowchart for Step 5



25

 **Unilateral Hip Rotation Range of Motion Asymmetry in Patients With Sacroiliac Joint Regional Pain**
Cibulka et al. *Spine*: Volume 23(9) 1 May 1998 pp 1009-15

- **The patients with LBP with evidence of SIJ dysfunction had significantly reduced hip internal rotation unilaterally, specifically on the side of the posterior innominate.**
- “The presence of such asymmetry in patients with low back pain may help identify those with sacroiliac joint dysfunction.”



26



“The Hip Joint: Myofascial and Joint Patterns,”
Marc Heller, DC, *Dynamic Chiropractic*, May 7, 2007:




- **Internal rotation is the first motion affected with hip capsular problems and DJD.**
- **“When you find a severe lack of internal rotation in the hip and this motion cannot be restored, the patient is usually on their way to hip replacement surgery.”**

ImpulseAdjusting
TECHNIQUE


27

Example: Reduced
internal rotation
of right hip:


Step 5




28



29




30



Posterior Hip Adjustment for IR Dysfunction


Posterior Hip Adjustment

SCP: Posterior aspect of the greater trochanter
LOD: Anterior
Setting: High or Medium




ImpulseAdjusting
TECHNIQUE


31




32



Step 5: Hip Rotation Tests




Internal Rotation of the Left Hip




External Rotation of the Left Hip

Note asymmetries



33

Example: Reduced external rotation of right hip:



Step 5

Hip Rotation Tests

50-60 degrees

External Rotation

Positive

Hyperactive or Inhibited Internal Rotators

Adjust Anterior Trochanter

34

Hip Internal Rotators: TFL
Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

ImpulseAdjusting
TECHNIQUE™

This slide features a header with a small image of a handheld device and a title 'Hip Internal Rotators: TFL' with a reference to Kendall et al (2005). Below the header, there is a large anatomical diagram of a right leg in a flexed position, with the Tensor Fasciae Latae (TFL) muscle highlighted in black. An inset photograph shows a person's leg in a similar flexed position, with hands placed on the thigh and knee to demonstrate the muscle's location.

35

Hip Internal Rotators: TFL
Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

Anterior fibers
Gluteal aponeurosis
GLUTEUS MEDIUS

GLUTEUS MINIMUS

ImpulseAdjusting
TECHNIQUE™


This slide features a header with a small image of a handheld device and a title 'Hip Internal Rotators: TFL' with a reference to Kendall et al (2005). Below the header, there are two anatomical diagrams of a right leg. The left diagram shows the Gluteus Medius muscle, with labels for 'Anterior fibers' and 'Gluteal aponeurosis'. The right diagram shows the Gluteus Minimus muscle. The slide also includes the 'ImpulseAdjusting TECHNIQUE' logo.

36

External Rotators
Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS



Posterior fibers

Gluteal aponeurosis

GLUTEUS MEDIUS

ImpulseAdjusting
TECHNIQUE™

37

Anterior Hip Adjustment
 (For Loss of Hip External Rotation)

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

Anterior Hip Adjustment

SCP: Anterior aspect of the greater trochanter

LOD: Posterior



Setting: High (or Medium in special circumstances)



ImpulseAdjusting
TECHNIQUE™

38

Step 6:
Lower Extremity Compression Tests






Compression Tests for Lateral Hamstring and Peroneus Longus:
Pressure applied along the muscle in a superior to inferior sweeping motion searching for active trigger points, identified by the “jump sign.”




39

Lateral Hamstring (Biceps Femoris m.)
Muscles, Testing and Function with Posture and Pain, Kendall et al. 2005, Lippincott Williams & Wilkins

Short head Long head Sacrotuberous lig.



40




Step 6. Lower Extremity Compression Test
(for Lateral Hamstring and Peroneal m. involvement)

mechanical INNOVATIONS


Impulse Adjusting TECHNIQUE

41




42

Step 6



43




Lateral Hamstring Adjustments



Lateral Hamstring Adjustments

SCP: Taught/Tender area of muscle belly of Biceps Femoris m.
LOD: Anterior
Setting: Medium or Low




ImpulseAdjusting
TECHNIQUE

44

Step 6




Peroneal Compression Test



45

Peroneal m.


Muscles, Testing and Function with Posture and Pain, Kendall et al, 2005, Lippincott Williams & Wilkins



PERONEUS LONGUS

PERONEUS BREVIS

46




Peroneal m. Adjustments

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS

Peroneal m. Adjustments

SCP: Taught/Tender area of muscle belly of Peroneal m.
LOD: Anterior
Setting: Medium or Low



ImpulseAdjusting
TECHNIQUE™

47



Step 7: Supra-Pelvic Compression Test

ImpulseAdjusting
TECHNIQUE™

mechanical
INNOVATIONS



Compression test on the quadratus lumborum

ImpulseAdjusting
TECHNIQUE™

48




Step 7. Supra-Pelvic Compression Tests
(for Quadratus Lumborum m. involvement)




ImpulseAdjusting
TECHNIQUE

49



50

Step 7

Supra-pelvic
Compression
Tests


Quadratus
Lumborum

Positive Findings:
Myospasm
Tenderness
Trigger Points
Guarding

Positive

**Hyperactive
Quadratus**

Impulse the
Quadratus
Lumborum



51

 Step 7. Supra-Pelvic Compression Tests
(for Quadratus Lumborum m. involvement)

ImpulseAdjusting
TECHNIQUE

mechanical
INNOVATIONS



ImpulseAdjusting
TECHNIQUE

52



Quadratus Lumborum m. Adjustment

mechanical
INNOVATORS

Quadratus Lumborum m. Adjustment

SCP: Taught/Tender area of muscle belly of Quadratus Lumborum m.

LOD: Anterior

Setting: Medium or Low



ImpulseAdjusting
TECHNIQUE

53

Adjustment Documentation Example

| Adjust | Area | Left | Right | Bilat... | Technique | Indicated by | Tende... | Modifier | Quality |
|-------------------------------------|-------------------|-------------------------------------|-------------------------------------|-------------------------------------|------------|--|-------------------------------------|----------|------------------|
| <input checked="" type="checkbox"/> | Occiput | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Impulse IQ | Indicated by reduced Cervical Rotation | <input type="checkbox"/> | | |
| <input type="checkbox"/> | C1 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | C2 | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Impulse IQ | Indicated by reduced Cervical Rotation | <input type="checkbox"/> | | |
| <input type="checkbox"/> | C3 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | C4 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Cervical Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | C5 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Cervical Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | C6 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Cervical Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | C7 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Cervical Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | T1 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Thoracic Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | T2 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Thoracic Spinous Compression | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T3 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T4 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T5 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T6 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | T7 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Thoracic Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | T8 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Thoracic Spinous Compression | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T9 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T10 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T11 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | T12 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | L1 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | L2 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | L3 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Lumbar Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | L4 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Lumbar Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | L5 | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by positive Lumbar Spinous Compression | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | Right SI | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Impulse IQ | Indicated by pain on palpation | <input checked="" type="checkbox"/> | mild | (quality: sharp) |
| <input checked="" type="checkbox"/> | Left SI | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | Impulse IQ | Indicated by pelvic asymmetry | <input checked="" type="checkbox"/> | mild | (quality: sharp) |
| <input checked="" type="checkbox"/> | Sacrum | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | Impulse IQ | Indicated by pelvic asymmetry | <input type="checkbox"/> | | |
| <input type="checkbox"/> | Coccyx | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input type="checkbox"/> | Scapula | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |
| <input checked="" type="checkbox"/> | Trochanter | <input type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | Impulse IQ | Indicated by joint hypomobility | <input type="checkbox"/> | | |
| <input type="checkbox"/> | Quadratus Lumb... | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | <input type="checkbox"/> | | |

54




The Lumbar Spine Region



NEUROMECHANICAL INNOVATIONS

55

Lumbar Zygapophyseal Joint Referral Patterns (Continued)



Trochanteric region

Posterior thigh region

Lum













Groin region

Figure 4-17
Lumbar zygapophyseal joint pain referral patterns. Zygapophyseal pain patterns of the lumbar spine as described by Fukui and colleagues.¹⁰ Lumbar zygapophyseal joints L1/2, L2/3, and L4/5 always referred pain to the lumbar spine region. Primary referral to the gluteal region was from L5/S1 (68% of the time). Levels L2/3, L3/4, L4/5, and L5/S1 occasionally referred pain to the trochanteric region (10% to 16% of the time). Primary referral to lateral thigh, posterior thigh, and groin regions were most often from L3/4, L4/5, and L5/S1 (5% to 30% of the time).

56

Neurological Examination

Diagnostic Utility of the Sensation, Manual Muscle Testing, and Reflexes for Lumbar Spinal Stenosis

| Level of Herniation | Pain | Numbness | Weakness | Atrophy | Reflexes |
|---|--|--|---|--|--|
|  <p>L3-4 disc; 4th lumbar nerve root</p> |  <p>Lower back, hip, posterolateral thigh, anterior leg</p> |  <p>Anteromedial thigh and knee</p> |  <p>Quadriceps</p> |  <p>Quadriceps</p> |  <p>Knee jerk diminished</p> |
|  <p>L4-5 disc; 5th lumbar nerve root</p> |  <p>Over sacroiliac joint, hip, lateral thigh, and leg</p> |  <p>Lateral leg, web of great toe</p> |  <p>Dorsiflexion of great toe and foot; difficulty walking on heels; foot drop may occur</p> |  <p>Minor</p> |  <p>Changes uncommon (absent or diminished posterior tibial reflex)</p> |

57










| | | | | | |
|--|--|---|--|--|--|
|  <p>L5-S1 disc; 1st sacral nerve root</p> |  <p>Over sacroiliac joint, hip, posterolateral thigh, and leg to heel</p> |  <p>Back of calf; lateral heel, foot and toe</p> |  <p>Plantar flexion of foot and great toe may be affected; difficulty walking on toes</p> |  <p>Gastrocnemius and soleus</p> |  <p>Ankle jerk diminished or absent</p> |
|  <p>Massive midline protrusion</p> | <p>Lower back, thighs, legs, and/or perineum depending on level of lesion; may be bilateral</p> | <p>Thighs, legs, feet, and/or perineum; variable; may be bilateral</p> | <p>Variable paralysis or paresis of legs and/or bowel and bladder incontinence</p> | <p>May be extensive</p> |  <p>Ankle jerk diminished or absent</p> |

Figure 4-20

Clinical features of herniated lumbar nucleus pulposus.

58



59

Postural Assessment

Reliability of Postural Assessment

| Test and Study | Description and Positive Findings | Population | Inter-examiner Reliability |
|---|---|--|-----------------------------|
| Forward head ¹⁴ | "Yes" if the patient's external auditory meatus was anteriorly deviated (anterior to the lumbar spine) | 22 patients with mechanical neck pain | $\kappa = -.10 (-.20, .00)$ |
| Excessive shoulder protraction ¹⁴ | "Yes" if the patient's acromions were anteriorly deviated (anterior to the lumbar spine) | | $\kappa = .83 (.51, 1.0)$ |
| C7-T2 excessive kyphosis ¹⁴ | Recorded as "normal" (no deviation), "excessive kyphosis," or "diminished kyphosis." Excessive kyphosis was defined as an increase in the convexity, and diminished kyphosis was defined as a flattening of the convexity of the thoracic spine (at each segmental group) | | $\kappa = .79 (.51, 1.0)$ |
| T3-5 excessive kyphosis ¹⁴ | | | $\kappa = .69 (.30, 1.0)$ |
| T3-5 decreased kyphosis ¹⁴ | | | $\kappa = .58 (.22, .95)$ |
| T6-10 excessive kyphosis ¹⁴ | | | $\kappa = .90 (.74, 1.0)$ |
| T6-10 decreased kyphosis ¹⁴ | | | $\kappa = .90 (.73, 1.0)$ |
| Kyphosis ¹¹ | With patient standing, examiner inspects posture from the side. Graded as "present" or "absent" | 111 adults age ≥ 60 with chronic low back pain and 20 asymptomatic patients | $\kappa = .21$ |
| Scoliosis ¹¹ | With patient standing, examiner runs finger along spinous processes. Patient bends over and examiner assesses height of paraspinal musculature. Graded as "present" or "absent" | | $\kappa = .33$ |
| Functional leg length discrepancy ¹¹ | Compare bilateral iliac crest height with patient standing. Graded as "symmetrical" or "asymmetrical" | | $\kappa = .00$ |

ICC or κ Interpretation
 .81-1.0 Substantial agreement
 .61-.80 Moderate agreement
 .41-.60 Fair agreement
 .11-.40 Slight agreement
 0-.10 No agreement

60

Passive Intervertebral Motion

Reliability of Assessing Limited or Excessive Passive Intervertebral Motion

| ICC or κ | Interpretation |
|-----------------|-----------------------|
| .81-1.0 | Substantial agreement |
| .61-.80 | Moderate agreement |
| .41-.60 | Fair agreement |
| .11-.40 | Slight agreement |
| .0-.10 | No agreement |

| Test and Study | Description and Positive Findings | Population | Reliability |
|--|---|--|---|
| Upper lumbar segmental mobility ³⁶ | With patient prone, examiner applies a posteroanterior force to the spinous process and lumbar facets of each lumbar vertebra. Mobility of each segment is judged as "normal" or "restricted" | 39 patients with low back pain | (Spinous) Inter-examiner $\kappa = .02$ (-.27, .32) (Left facet) Inter-examiner $\kappa = .17$ (-.14, .48) (Right facet) Inter-examiner $\kappa = -.01$ (-.33, .30) |
| Lower lumbar segmental mobility ³⁶ | | | (Spinous) Inter-examiner $\kappa = -.05$ (-.36, .27) (Left facet) Inter-examiner $\kappa = -.17$ (-.41, .06) (Right facet) Inter-examiner $\kappa = -.12$ (-.41, .18) |
| Identifying the least mobile segment ³⁷ | With patient prone, examiner applies a posteroanterior force to the spinous process of each lumbar vertebra | 29 patients with central low back pain | Inter-examiner $\kappa = .71$ (.48, .94) |
| Identifying the most mobile segment ³⁷ | | | Inter-examiner $\kappa = .29$ (-.13, .71) |
| Posterior to anterior (PA) stiffness ³⁸ | | | Intra-examiner $\kappa = .54$ Intra-examiner (± 1 level) $\kappa = .64$ Inter-examiner $\kappa = .23$ Inter-examiner (± 1 level) $\kappa = .52$ |

61

Passive Intervertebral Motion (continued)

Reliability of Assessing Limited or Excessive Passive Intervertebral Motion

| ICC or κ | Interpretation |
|-----------------|-----------------------|
| .81-1.0 | Substantial agreement |
| .61-.80 | Moderate agreement |
| .41-.60 | Fair agreement |
| .11-.40 | Slight agreement |
| .0-.10 | No agreement |

| Test and Study | Description and Positive Findings | Population | Reliability |
|---|---|---|---|
| Determination of segmental fixations ⁴⁰ | Passive motion palpation is performed and the segment is considered fixated if a hard end-feel is noted during the assessment | 60 asymptomatic volunteers | Intra-examiner κ ranged from -.09 to .39 Inter-examiner κ ranged from -.06 to .17 |
| Passive motion palpation ⁴¹ | | 21 symptomatic and 25 asymptomatic subjects | Inter-examiner $\kappa =$ ranged from -.03 to .23 with a mean of .07 |
| Segmental mobility testing ⁴² | With patient side-lying with hips and knees flexed, examiner assesses mobility while passively moving the patient. Examiner determines whether mobility of the segment is "decreased," "normal," or "increased" | 71 patients with low back pain | Inter-examiner $\kappa = .54$ |
| Hypermobility at any level ⁴³ | With patient prone, examiner applies a posteroanterior force to the spinous process of each lumbar vertebra. Mobility of each segment is judged as "normal," "hypermobile," or "hypomobile" | 49 patients with low back pain referred for flexion-extension radiographs | Inter-examiner $\kappa = .48$ (.35, .61) |
| Hypomobility at any level ⁴³ | | | Inter-examiner $\kappa = .38$ (.22, .54) |
| Determination of posteroanterior spinal stiffness ⁴³ | Five raters tested lumbar spinal levels for posteroanterior mobility and graded each on an 11-point scale ranging from "markedly reduced stiffness" to "markedly increased stiffness" | 40 asymptomatic individuals | Inter-examiner ICC in the first study = .55 (.32, .79) Inter-examiner ICC in the second study = .77 (.57, .89) |
| Posteroanterior mobility testing ⁴⁴ | With the patient prone, examiner evaluates posteroanterior motion mobility. Mobility is scored on a 9-point scale ranging from "severe excess motion" to "no motion," and the presence of pain is recorded | 18 patients with low back pain | Inter-examiner ICC = .25 (.00, .39) |


62

| | | | |
|---|--|---|--|
| Posteroanterior mobility testing ⁴⁴ | With the patient prone, examiner evaluates posteroanterior motion mobility. Mobility is scored on a 9-point scale ranging from "severe excess motion" to "no motion," and the presence of pain is recorded | 18 patients with low back pain | Inter-examiner ICC = .25 (.00, .39) |
| Segmental mobility testing ⁴⁵ | With patient prone, examiner applies an anteriorly directed force over the spinous process of the segment to be tested. Examiner grades the mobility as "hypermobile," "normal" or "hypomobile" | 63 patients with current low back pain | Inter-examiner κ ranged from -.20 to .26 depending on level tested |
| Identification of a misaligned vertebra ⁴¹ | Static palpation is used to determine the relationship of one vertebra to the vertebra below | 21 symptomatic and 25 asymptomatic subjects | Inter-examiner κ ranged from -.04 to .03 with a mean of .00 |
| Detection of a segmental lesion T11-L5/S1 ⁴⁶ | Two clinicians used visual postural analysis, pain descriptions, leg length discrepancy, neurological examination, motion palpation, static palpation, and any special orthopaedic tests to determine the level of segmental lesion. | 19 patients with chronic mechanical low back pain | Intra-examiner κ = -.08 to .43 Inter-examiner κ = -.16 to .25 |

63

Passive Intervertebral Motion

Reliability of Assessing Painful Passive Intervertebral Motion



| ICC or κ | Interpretation |
|-----------------|-----------------------|
| .81-1.0 | Substantial agreement |
| .61-.80 | Moderate agreement |
| .41-.60 | Fair agreement |
| .11-.40 | Slight agreement |
| .0-.10 | No agreement |

Figure 4-26
Assessment of posteroanterior segmental mobility.

64

Figure 4-26

Assessment of posteroanterior segmental mobility.

| Test and Study | Description and Positive Findings | Population | Reliability | |
|---|--|--|---|---------------------------|
| | | | Intra-examiner | Inter-examiner |
| Spring test T10-T7 ⁴⁷ | With patients in the prone position the therapist applies a posteroanterior force to the spinous processes of T7-L5. The pressure of each force is held for 20 seconds. Considered positive if the force produces pain | 84 subjects, of whom 53% reported experiencing low back symptoms within the last 12 months | $\kappa = .73 (.39-.1.0)$ | $\kappa = .12 (-.18-.41)$ |
| Spring test L2-T11 ⁴⁷ | | | $\kappa = .78 (.49-.1.0)$ | $\kappa = .36 (.07-.66)$ |
| Spring test L5-L3 ⁴⁷ | | | $\kappa = .56 (.18-.94)$ | $\kappa = .41 (.12-.70)$ |
| Pain with upper lumbar mobility testing ²⁸ | With patient prone, examiner applies a posteroanterior force to the spinous processes and lumbar facets of each lumbar vertebra. Response at each segment is judged as "painful" or "not painful" | 39 patients with low back pain | (Spinous) Inter-examiner $\kappa = .21 (-.10-.53)$ (Left facet) Inter-examiner $\kappa = .46 (.17-.75)$ (Right facet) Inter-examiner $\kappa = .38 (.06-.69)$ | |
| Pain with lower lumbar mobility testing ²⁸ | | | (Spinous) Inter-examiner $\kappa = .57 (.32-.83)$ (Left facet) Inter-examiner $\kappa = .73 (.51-.95)$ (Right facet) Inter-examiner $\kappa = .52 (.25-.79)$ | |
| Pain provocation ⁴⁵ | With patient prone, examiner applies an anteriorly directed force over the spinous processes of the segment to be tested. Considered positive if pain is reproduced | 63 patients with current low back pain | Inter-examiner κ ranged from .25 to .55 depending on the segmental level tested | |
| Pain during mobility testing ³¹ | | 49 patients with low back pain referred for flexion-extension radiographs | Inter-examiner $\kappa = .57 (.43-.71)$ | |

LITERATURE REVIEW

MANUAL EXAMINATION OF THE SPINE: A SYSTEMATIC CRITICAL LITERATURE REVIEW OF REPRODUCIBILITY

Mette Jensen Stochkendahl, DC,^a Henrik Wulff Christensen, DC, MD, PhD,^b Jan Hartvigsen, DC, PhD,^c Werner Vach, PhD,^d Mitchell Haas, DC, MA,^e Lise Hestbaek, DC, PhD,^f Alan Adams, DC, MS, MEd,^g and Gert Bronfort, DC, PhD^h


ABSTRACT

Objective: Poor reproducibility of spinal palpation has been reported in previously published literature, and authors of recent reviews have posted criticism on study quality. This article critically analyzes the literature pertaining to the inter- and intraobserver reproducibility of spinal palpation to investigate the consistency of study results and assess the level of evidence for reproducibility.


Methods: Systematic review and meta-analysis were performed on relevant literature published from 1965 to 2005, identified using the electronic databases MEDLINE, MANTIS, and CINAHL and checking of reference lists. Descriptive data from included articles were extracted independently by 2 reviewers. A 6-point scale was constructed to assess the methodological quality of original studies. A meta-analysis was conducted among the high-quality studies to investigate the consistency of data, separately on motion palpation, static palpation, osseous pain, soft tissue pain, soft tissue changes, and global assessment. A standardized method was used to determine the level of evidence.

Results: The quality score of 48 included studies ranged from 0% to 100%. There was strong evidence that the interobserver reproducibility of osseous and soft tissue pain is clinically acceptable ($\kappa \geq 0.4$) and that intraobserver reproducibility of soft tissue pain and global assessment are clinically acceptable. Other spinal procedures are either not reproducible or the evidence is conflicting or preliminary. (*J Manipulative Physiol Ther* 2006;29:475-485)

Key Indexing Terms: Reproducibility of Results; Palpation; Literature Review; Diagnostic Tests; Spine; Meta-Analysis




67




68




69




70




71



72



73



74



Lower Thoracic Spinous Compression Test:

Impulse Adjusting
TECHNIQUE

mechanical
INNOVATIONS

- Proceed from lumbar spine up the thoracic spine to the apex of the thoracic curve (usually around T6-7)
- Rationale for spinous compression (producing posterior to anterior movement) is that the zygapophyseal joints in the lower thoracic spine are in the sagittal plane (like the lumbar spine facets) allowing mostly flexion.
- Amount of flexion per segment increases from approx 4 degrees at T6-7 to 12 degrees at T11-12.
- Rotation is very limited, especially in the lowest levels.

Impulse Adjusting
TECHNIQUE